

Analog EMS Unit

ReliaMed® products offer a wide variety of choices, and are manufactured with a focus on quality — at a value that's hard to ignore.



ZZA550
HCPCS E0745

ReliaMed® Analog EMS Unit

- Features protective cap over the amplitude adjustment knobs to prevent unintended changes
- Comes in a hard-sided, protective, plastic case
- Sold by the each or 30 per case
- Utilizes standard lead wires and "pigtail"-style electrodes
- Includes two lead wires, one pack (four) electrodes, manual and 9-volt battery

ZZA550 Analog EMS Unit (HCPCS E0745)	
Channel	Dual, isolated between channels
Modes	Synchronous, constant and alternation
Pulse Intensity	Adjustable 0 – 80 mA peak into 500 Ω load each channel, constant current
Wave Form	Asymmetrical, biphasic, square pulse
Pulse Rate	2 – 120 Hz (adjustable)
Pulse Width	250 μ s (fixed)
ON Ramp	1 – 8 seconds (adjustable)
Cycle ON and OFF Times	1 – 30 seconds (adjustable)
Timer	30- or 60-minute; or continuous

Call your **ReliaMed** Sales Representative for more information.

Toll-free **1-800-409-2848** Fax **330-405-7117** Web **www.reliamedproducts.com**