

Compressor Nebulizer

ReliaMed® products offer a wide variety of choices, and are manufactured with a focus on quality — at a value that's hard to ignore.

ReliaMed® Compressor Nebulizer

- For the successful treatment of asthma, allergies and other respiratory disorders
- Creates a stream of air that travels through the clear tube to the nebulizer, which converts medication into an aerosol mist for easy inhalation
- 5-mL medication capacity
- Operable up to a 45° angle
- Compact and lightweight
- Easy-to-use, one-button operation
- Only uses ReliaMed® Air Filters
- Ideal for all ages
- Includes piston-style compressor, nebulizer, 7' tubing, angled mouthpiece, five filters, and detailed guidebook in English and Spanish
- Three-year limited warranty



| Item# | Description | Sold By | Case Qty. | HCPCS |
|---------|-----------------------|---------|-----------|-------|
| ZRCN01 | Compressor Nebulizer | each | 6 | E0570 |
| ZRN01RF | Nebulizer Air Filters | 10/pack | 100 | A7013 |

Call your **ReliaMed** Sales Representative for more information.

Toll-free **1-800-409-2848** Fax **330-405-7117** Web **www.reliamedproducts.com**