



Digital Moist Heating Pads

ReliaMed® products offer a wide variety of choices, and are manufactured with a focus on quality — at a value that’s hard to ignore.

ReliaMed® Digital Moist Heating Pads

- Controlled moist heat for temporary pain relief associated with muscle strains and sprains, joint pain and inflammation
- Temporarily ease back pain caused by overexertion, stress and strain
- Do not require water – produce their own moisture using humidity in the air
- Easy-to-read LCD controls allow user to set and monitor treatment temperature and time
- Lock-out feature guards against unintended touching of the controls during treatment
- Shut-off can be set in one-minute increments up to one hour
- Automatic shut-off after 15 minutes
- Include heating pad, fleece cover and instruction manual
- One-year warranty



Item#	Size	Sold By	Case Qty.	HCPCS
ZHP147	14" x 7"	each	12	E0215
ZHP1414	14" x 14"	each	12	E0215
ZHP1427	14" x 27"	each	12	E0215
ZHP832U	8" x 32", U-Shaped	each	12	E0215

Call your **ReliaMed** Sales Representative for more information.

Toll-free **1-800-409-2848** Fax **330-405-7117** Web **www.reliamedproducts.com**