

Insertion Trays

ReliaMed® products offer a wide variety of choices, and are manufactured with a focus on quality — at a value that's hard to ignore.

ReliaMed® Foley Catheter Insertion Trays

- Contain all of the components required to insert and flush a foley catheter, according to physician's instruction
- Include 1,000-mL tray, prefilled sterile water syringe, three PVP swabs, underpad, fenestrated drape, 5 cc of lubricating jelly, one pair of PVC gloves, specimen container and label
- Sterile

latex free

Item#	Description	Sold By	Case Qty.	HCPCS
ZRCIT10CC	With 10-cc Prefilled Syringe	each	20	A4310
ZRCIT30CC	With 30-cc Prefilled Syringe	each	20	A4310



ReliaMed® Foley Catheter Insertion Trays Compare To:	
Item#	Vendor & Product
57782100	Bard® Universal Catheter Insertion Tray, 10-cc Syringe
57782102	Bard® Universal Catheter Insertion Tray, 30-cc Syringe
6876010	Covidien KENGUARD™ Universal Catheter Insertion Tray, 10-cc Syringe
6876030	Covidien KENGUARD™ Universal Catheter Insertion Tray, 30-cc Syringe
RU76710	Rüsch Foley Catheter Procedure Tray, 10-cc Syringe
RU76700	Rüsch Foley Catheter Procedure Tray, 30-cc Syringe

Bard® is a registered trademark of C.R. Bard, Inc.;
KENGUARD™ is a trademark of Covidien AG

Call your **ReliaMed** Sales Representative for more information.
Toll-free **1-800-409-2848** Fax **330-405-7117** Web **www.reliamedproducts.com**