


Red Rubber Catheters

ReliaMed® products offer a wide variety of choices, and are manufactured with a focus on quality — at a value that's hard to ignore.

ReliaMed® Red Rubber Catheters

- Feature smooth, rounded tips for ease of insertion
- Two opposing, smooth-cut eyes allow for high drainage flow, and to aspirate clots and debris from the bladder
- Constructed of heavy-duty latex for additional stiffness, which aids in insertion
- For single use only; sterile
- 16"
- Contain latex

Item#	Size	Sold By	HCPCS	Case Qty.
ZRRRC8	8 fr	each	A4351	200
ZRRRC10	10 fr	each	A4351	200
ZRRRC12	12 fr	each	A4351	200
ZRRRC14	14 fr	each	A4351	200
ZRRRC16	16 fr	each	A4351	200
ZRRRC18	18 fr	each	A4351	200
ZRRRC20	20 fr	each	A4351	200
ZRRRC22	22 fr	each	A4351	200



ReliaMed® Red Rubber Catheters Compare To:	
Item#	Vendor & Product
59408-22	Bard® Red Rubber Catheters, 8 fr – 22 fr (even numbers)
RU35108-22	Rüsch Red Rubber Robinson/Nelaton Catheters, 8 fr – 22 fr (even numbers)
61660143	Covidien Dover™ Red Rubber Robinson Catheters, 14 fr
61660168	Covidien Dover™ Red Rubber Robinson Catheters, 16 fr
61660218	Covidien Dover™ Red Rubber Robinson Catheters, 20 fr

Bard® is a registered trademark of C.R. Bard, Inc.;
Dover™ is a registered trademark of Covidien AG

Call your **ReliaMed** Sales Representative for more information.
Toll-free **1-800-409-2848** Fax **330-405-7117** Web **www.reliamedproducts.com**