

Portable Ultrasound

ReliaMed® products offer a wide variety of choices, and are manufactured with a focus on quality — at a value that's hard to ignore.

ReliaMed® Portable Ultrasound

- Portable ultrasound unit features three power settings: low, medium and high
- Generates deep heat within body tissues for treatment of selected medical conditions, such as pain relief, muscle spasms and joint contractures
- Sold by the each or 20 per case
- Comes in a hard-sided, protective, plastic case
- Includes AC adapter, ultrasound gel, manual and carrying case
- Caution: This product should not be used for the treatment of malignancies



ZZA1000
HCPCS E1399
Portable Ultrasound

ZZA1000 Portable Ultrasound (HCPCS E1399)	
Waveform	Pulse
Output Frequency	1 MHz ± 10%
Pulse Width	2 ms ± 10%
Repetition Rate	150 Hz ± 10%
Output Intensities	Low, medium and high
Temporal Maximum Power	4 W ± 10%
Timer	30-minute with automatic shut-off

Call your **ReliaMed** Sales Representative for more information.

Toll-free **1-800-409-2848** Fax **330-405-7117** Web **www.reliamedproducts.com**